

LF MOBILE X-RAY : 100mA

X-Ray Power, Unmatched Clarity



*Applied For / Not For Sale


Technical Highlights

POWER RATING (KW)	8KW
POWER SUPPLY	Single Phase 230V, 50Hz, Line Resistance 0.40 max
RADIOLOGY	100mA at 100kVp
X-RAY GENERATOR	Single Phase (2 pulse), Full wave rectified
TIMER (S)	0.06-5
TUBE CURRENT (MA)	25 : 50 : 75 : 100
TUBE VOLTAGE (KVP)	45-100 kVp (15 steps)

Clinical Images*



NOTE: *For reference purpose only



Edusoft Healthcare Pvt. Ltd.

B-83, Phase-II, Mangolpuri Industrial Area, New Delhi-110034, India

ceo@edusofthealth.com

1800-120-280-280

+91 8851248073, +91 8826066621

- Noida • Lucknow • Chandigarh • Pune • Bhopal • Chennai • Bengaluru • Kolkata • Patna • Bhopal
- Hyderabad • Jaipur • Guwahati • Ahmedabad • Mumbai

Copyright © 2024 Edusoft Healthcare Pvt. Ltd. All rights reserved. Specifications and designs are subject to change without notice. Non-metric weights and measurements are approximate. All data were deemed correct at time of creation. Edusoft is not liable for errors or omissions. All brand, product, service names and logos are trademarks and/or registered trademarks of their respective owners and are hereby recognized and acknowledged.



Scan Code or visit
www.edusofthealth.com
 to learn more