

ERAY™

ULTRA PORTABLE HANDHELD X-RAY

Ultra-Light, Powerful & Affordable



Mass Testing
For TB



Weight
Under 2kg



Lithium
Battery Powered



LED Collimator



150 Exposures
per Charge



C-Type
Charger



Technical Highlights

MODEL	ERAY SMART 3HS
APPLICATION	Chest X-Ray + Extremities
WORKING FREQUENCY	50kHz
NOMINAL POWER	210W
kV RANGE	40kV ~ 70kV
mA RANGE	1mA ~ 3mA (1 Step)
EXPOSURE TIME	40ms ~ 3200ms
FOCAL SPOT	0.4mm
HEAT STORAGE CAPACITY	4.5kJ
WEIGHT	Less than 2kg

X-Ray Solution Kit Available



ERAY SMART 3HS


TRIPOD STAND
(Optional)


Clinical Images*





*For reference purpose only

Edusoft Healthcare Pvt. Ltd.

 B-83, Phase-II, Mangolpuri Industrial Area, New Delhi-110034, India

 ceo@edusofthealth.com

 1800-120-280-280

 +91 8851248073, +91 8826066621

• Noida • Lucknow • Chandigarh • Pune • Bhopal • Chennai • Bengaluru • Kolkata
• Hyderabad • Jaipur • Guwahati • Ahmedabad • Mumbai

Copyright © 2025 Edusoft Healthcare Pvt. Ltd. All rights reserved. Specifications and designs are subject to change without notice. Non-metric weights and measurements are approximate. All data were deemed correct at time of creation. Edusoft is not liable for errors or omissions. All brand, product, service names and logos are trademarks and/or registered trademarks of their respective owners and are hereby recognized and acknowledged.



Scan Code or visit
www.edusofthealth.com
to learn more