



## *Patient Care Furniture*



ISO CERTIFIED 13485:2016 & 9001:2015



Indtech expertise in manufacturing innovative, comfortable and easy-to-use healthcare furnitures with utmost safety for patients in mind. Working for frontline healthcare workers including doctors and nurses and patients who seeks care and comfort in their difficult times, Indtech never compromise quality.

Our products ranges from ICU advanced care and ICU support to Dialysis, Examination support, Ward care, General support, Transfer support and Mother care. The right talent in production applied upfront allows us to meet almost all of our client's requirements in time.

## QUALITY POLICY

The quality policy of IndTech Health Care Pvt. Ltd., is to provide products that meet or exceed customer expectations, comply with all applicable regulations, continually improve in everything we do, and maintain the effectiveness of the quality management system.

### WHAT MAKES US SPECIAL



#### CUSTOMISED PRODUCTS

We provide our customers with the products that meet their specific demands and requirements.



#### ISO 13485:2016 CERTIFIED

The products are accredited to ISO 13485:2016 & 9001: 2015 and CE certification for European Economic Area



#### INTERNATIONAL STANDARDS

Every time we make sure that the products we manufacture meets the international export standards.



#### PRODUCT MAINTENANCE

We don't quit once the product is handed over to our clients. We provide maintenance services for the products on site.



#### ON TIME DELIVERY

We make sure that the products we manufacture get delivered to our clients on time.



#### STATE-OF-THE-ART PRODUCTION SYSTEM

We use the latest technologies to manufacture the products and ensure the utmost quality.

In the face of laborious working environment in hospitals, Indtech offers comfortable and easy to handle furnitures with utmost safety for patients. With varieties of different designs and grades from deluxe to normal, we meet requirements for each and every hospitals. Check out our products to know more.





# LEARN MORE ABOUT INDTECH PATIENT CARE FURNITURES.

---


ICU Advanced Care  
Dialysis  
Ward Care  
Transfer Support  
Mother Care

ICU Support  
Examination Support  
General Support



### IHC 1201 – TRANSFER TROLLEY SYSTEM

- **Back-rest elevation** controlled by ratchet.
- **IV pole:** Provision at two diagonal ends and made of stainless steel.
- **Removable stretcher** made of stainless steel.
- Two under carriages and one rolling stretcher.
- **Interlocking system** for undercarriages to roll over stretchers from one carriage to another.
- **Safety belts** provided.
- Pre-treated and epoxy powder coated.
- **Wheels:** Four 150mm non marking castors with diagonally opposite individual locking system.



### MATTRESS PLATFORM DIMENSIONS


Length : 2030mm  
Width : 565mm  
Height : 810mm

### ADJUSTABLE ANGLE RANGE (DEGREE)

Backrest angle : 0 - 85

### OPTIONAL ACCESSORIES

Mattress  
Telescopic SS Saline rod



Backrest



### IHC 1204 – EMERGENCY CASUALTY TROLLEY

- **Adjustable back-rest** controlled by ratchet.
- **Gas spring assisted lever** for controlling trendelenburg and reverse trendelenburg positions.
- **Bed top:** Made of radio permeable material for taking X-Rays.
- **X-Ray cassette holder** of full length under the bed top.
- **Side rails:** Swing-away type made of stainless steel.
- **IV rod:** made of stainless steel.
- **Oxygen cylinder provision** on the base.
- **Mechanical foot pedals** for hi-low positions assisted by hydraulic pump.
- Pre-treated and epoxy powder coated.
- **Wheels:** Four 150mm non marking castors with diagonally opposite individual locking system.

### MATTRESS PLATFORM DIMENSIONS

Length : 1830 mm

Width : 560mm

### HEIGHT ADJUSTMENT

Min : 750mm without mattress.


Max : 870mm without mattress.

### ADJUSTABLE ANGLE RANGE (DEGREE)

Backrest angle : 0 - 80

Trendelenburg angle : 0 - 20


Reverse Trendelenburg angle : 0 - 20




Backrest



Hi-Low



Trendelenburg



R.Trendelenburg



### IHC 1206 – EMERGENCY TROLLEY

- **Adjustable height** (scissor type).
- **Fixed stretcher** made of stainless steel.
- **Mounted on** four castors of 150mm (two with brakes).
- **Finish pre-treated** and epoxy powder coated.
- **Ratchet** backrest elevation.
- Telescopic ss saline rod.
- **Oxygen cylinder** provision.
- Stainless steel swing away side railings.
- Four rubber buffers provided to prevent damage to wall.
- **Foldable crank** system for hi-low operation.

### MATTRESS PLATFORM DIMENSIONS

Length : 1830 mm

Width : 560mm


### HEIGHT ADJUSTMENT

Min : 600mm without mattress.

Max : 810mm without mattress.

### ADJUSTABLE ANGLE RANGE (DEGREE)

Backrest angle : 0 - 80



Backrest



Hi-Low





#### **IHC 1207 – STRETCHER TROLLEY**

- Removable stretcher.
- Carriage mounted on four castors of 150mm (two with brakes).
- Finish pre-treated and epoxy powder coated.
- Oxygen cylinder provision.
- MS swing away side railings.
- IV pole provision at two diagonal ends.

#### **MATTRESS PLATFORM DIMENSIONS**

Length : 1830 mm  
Width : 560 mm  
Height : 810 mm



#### **IHC 1208 – STRETCHER TROLLEY (FULL SS)**

- Removable stretcher made of stainless steel.
- SS carriage mounted on four castors of 150mm (Two with brakes).
- IV pole provision at two diagonal ends.
- Telescopic SS saline rod.
- Oxygen cylinder provision.
- Stainless steel swing away side railings.
- Safety belts provided.

#### **MATTRESS PLATFORM DIMENSIONS**

Length : 1830 mm  
Width : 560 mm  
Height : 600 mm



#### **IHC 1209 – STRETCHER TROLLEY**

- Removable stretcher made of stainless steel.
- MS carriage mounted on four castors of 150mm (two with brakes).
- IV pole provision at two diagonal ends.
- Telescopic ss saline rod.
- Oxygen cylinder provision.
- MS swing away side railings.
- Safety belts provided.
- Finish pre-treated and epoxy powder coated.

#### **MATTRESS PLATFORM DIMENSIONS**

Length : 1830 mm  
Width : 560 mm  
Height : 600 mm



#### **IHC 1210 - STRETCHER TROLLEY MS**

- Removable stretcher made of MS.
- MS carriage mounted on four castors of 150mm (two with brakes).
- IV pole provision at two diagonal ends.
- Telescopic ss saline rod.
- Oxygen cylinder provision.
- MS swing away side railings.
- Safety belts provided.
- Finish pre-treated and epoxy powder coated.

#### **MATTRESS PLATFORM DIMENSIONS**

Length : 1830 mm  
Width : 560 mm  
Height : 600 mm



#### **IHC 1211 - BABY TROLLEY**

- Stretcher made of stainless steel
- Complete SS carriage mounted on four Castors of 100mm (Two with brakes)
- File tray
- Stainless steel swing away side railings
- Safety belts provided



#### **IHC 1212 - BABY TROLLEY**

- Frame work made of CRCA steel tube
- Mounted on four twin-wheel of 100mm.
- Finish pre-treated and epoxy powder coated.
- Baby mattress



#### **IHC1220 - SS WHEEL CHAIR**

- Framework made of stainless steel.
- Sitting position and back rest made of stainless steel mat finish sheet.
- Plastic folding foot step
- Two 150mm swivel castors at front
- Two solid mag wheels at back

#### **DIMENSIONS**

Length : 1120 mm

Width : 680 mm

Height : 935 mm



#### CORPORATE OFFICE

Indtech Healthcare Pvt. Ltd, No. 6/858-M,  
2nd Floor, Suite No.893, Valamkottil Towers  
Judgemukku, Kakkanad, Kochi - 682021  
Ph: +91 7034 111 182



#### CALICUT OFFICE

Indtech Healthcare Pvt. Ltd, Office 1513,  
5th Floor, HiLITE Business Park, Calicut- 673014  
Ph: +91 7560 822 822



#### FACTORY

Plot No: B5, Industrial Development Plot,  
Paloor Kotta, Palachode (P.O), Angadippuram (Via),  
Malappuram-679 338, Kerala, Ph: +91 9495 351458



ISO CERTIFIED 13485:2016 & 9001:2015



#### Customer Care

+91 7594 822 822

[www.indtechhc.com](http://www.indtechhc.com)

[info@indtechhc.com](mailto:info@indtechhc.com)