



Patient Care Furniture



ISO CERTIFIED 13485:2016 & 9001:2015



Indtech expertise in manufacturing innovative, comfortable and easy-to-use healthcare furnitures with utmost safety for patients in mind. Working for frontline healthcare workers including doctors and nurses and patients who seeks care and comfort in their difficult times, Indtech never compromise quality.

Our products ranges from ICU advanced care and ICU support to Dialysis, Examination support, Ward care, General support, Transfer support and Mother care. The right talent in production applied upfront allows us to meet almost all of our client's requirements in time.

QUALITY POLICY

The quality policy of IndTech Health Care Pvt. Ltd., is to provide products that meet or exceed customer expectations, comply with all applicable regulations, continually improve in everything we do, and maintain the effectiveness of the quality management system.

WHAT MAKES US SPECIAL



CUSTOMISED PRODUCTS

We provide our customers with the products that meet their specific demands and requirements.



ISO 13485:2016 CERTIFIED

The products are accredited to ISO 13485:2016 & 9001: 2015 and CE certification for European Economic Area



INTERNATIONAL STANDARDS

Every time we make sure that the products we manufacture meets the international export standards.



PRODUCT MAINTENANCE

We don't quit once the product is handed over to our clients. We provide maintenance services for the products on site.



ON TIME DELIVERY

We make sure that the products we manufacture get delivered to our clients on time.



STATE-OF-THE-ART PRODUCTION SYSTEM

We use the latest technologies to manufacture the products and ensure the utmost quality.

In the face of laborious working environment in hospitals, Indtech offers comfortable and easy to handle furnitures with utmost safety for patients. With varieties of different designs and grades from deluxe to normal, we meet requirements for each and every hospitals. Check out our products to know more.






LEARN MORE ABOUT INDTECH PATIENT CARE FURNITURES.

ICU Advanced Care
Dialysis
Ward Care
Transfer Support
Mother Care

ICU Support
Examination Support
General Support




IHC1401 – OBSTETRIC LABOR TABLE (FULLY SS)

- Headboard and hand grips made of stainless steel.
- Bed top: Three sections made of stainless steel.
- Lithotomy crutches with radial clamp.
- “U” cut on bed top for lithotomy and a stainless steel bowl under it.
- Leg section can be telescoped under the main frame for lithotomy.
- Adjustable back-rest controlled by ratchet.
- Crank system for trendelenburg and reverse trendelenburg positions.


MATTRESS PLATFORM DIMENSIONS

Length : 1890 mm
Width : 765 mm
Height : 800 mm




ADJUSTABLE ANGLE RANGE (DEGREE)


Backrest angle : 0 - 70
Trendelenburg angle : 0 - 20
Reverse Trendelenburg angle : 0 - 20



Backrest



Trendelenburg



R.Trendelenburg



IHC 1403 – LABOUR TABLE (FULLY SS)

- Bed top made of stainless steel
- Lithotomy crutches with radial clamp.
- “U” cut at the end for lithotomy.
- IV rod made of stainless steel.
- Available in 304 and 202 grades.

DIMENSIONS

Length : 1830 mm
Width : 765 mm
Height : 800 mm



IHC 1404 – LABOUR TABLE (FULLY SS)

- Bed top made of stainless steel
- Stainless steel lithotomy rods with Rexine angle straps.
- “U” cut at the end for lithotomy.
- IV rod made of stainless steel.
- Available in 304 and 202 grades.

DIMENSIONS

Length : 1830 mm
Width : 765 mm
Height : 800 mm



CORPORATE OFFICE

Indtech Healthcare Pvt. Ltd, No. 6/858-M,
2nd Floor, Suite No.893, Valamkottil Towers
Judgemukku, Kakkanad, Kochi - 682021
Ph: +91 7034 111 182



CALICUT OFFICE

Indtech Healthcare Pvt. Ltd, Office 1513,
5th Floor, HiLITE Business Park, Calicut- 673014
Ph: +91 7560 822 822



FACTORY

Plot No: B5, Industrial Development Plot,
Paloor Kotta, Palachode (P.O), Angadippuram (Via),
Malappuram-679 338, Kerala, Ph: +91 9495 351458



ISO CERTIFIED 13485:2016 & 9001:2015



Customer Care

+91 7594 822 822

www.indtechhc.com

info@indtechhc.com