



Patient Care Furniture



ISO CERTIFIED 13485:2016 & 9001:2015



Indtech expertise in manufacturing innovative, comfortable and easy-to-use healthcare furnitures with utmost safety for patients in mind. Working for frontline healthcare workers including doctors and nurses and patients who seeks care and comfort in their difficult times, Indtech never compromise quality.

Our products ranges from ICU advanced care and ICU support to Dialysis, Examination support, Ward care, General support, Transfer support and Mother care. The right talent in production applied upfront allows us to meet almost all of our client's requirements in time.

QUALITY POLICY

The quality policy of IndTech Health Care Pvt. Ltd., is to provide products that meet or exceed customer expectations, comply with all applicable regulations, continually improve in everything we do, and maintain the effectiveness of the quality management system.

WHAT MAKES US SPECIAL



CUSTOMISED PRODUCTS

We provide our customers with the products that meet their specific demands and requirements.



ISO 13485:2016 CERTIFIED

The products are accredited to ISO 13485:2016 & 9001: 2015 and CE certification for European Economic Area



INTERNATIONAL STANDARDS

Every time we make sure that the products we manufacture meets the international export standards.



PRODUCT MAINTENANCE

We don't quit once the product is handed over to our clients. We provide maintenance services for the products on site.



ON TIME DELIVERY

We make sure that the products we manufacture get delivered to our clients on time.



STATE-OF-THE-ART PRODUCTION SYSTEM

We use the latest technologies to manufacture the products and ensure the utmost quality.

In the face of laborious working environment in hospitals, Indtech offers comfortable and easy to handle furnitures with utmost safety for patients. With varieties of different designs and grades from deluxe to normal, we meet requirements for each and every hospitals. Check out our products to know more.





LEARN MORE ABOUT INDTECH PATIENT CARE FURNITURES.

ICU Advanced Care
Dialysis
Ward Care
Transfer Support
Mother Care

ICU Support
Examination Support
General Support



IHC 1700 – MOTORIZED EXAMINATION TABLE

- 75mm thick high density foam upholstered top in two sections.
- Adjustable Height and backrest on electric actuators.
- Provided with a drawers which can be opened from both sides of the couch for easy accessibility and paper roll holder at foot end of the couch.
- Chemically pre-treated and polyester epoxy powder coating finish.

MATTRESS PLATFORM DIMENSIONS

Length : 1830 mm

Width : 610 mm


HEIGHT ADJUSTMENT

Min : 485mm without mattress.

Max : 855mm without mattress.

ADJUSTABLE ANGLE RANGE (DEGREE)

Backrest angle : 0 - 80



Backrest



Hi-Low



IHC 1701 – EXAMINATION COUCH

- CRCA rectangular tube frame with CRCA sheet cabinets
- Retractable step stool
- Headrest adjustable on gas spring.
- Three cabinets and three drawers.
- Finish pre-treated and epoxy powder coated.
- B.P apparatus tray.

DIMENSIONS

Length : 1830 mm
Width : 610 mm
Height : 810 mm



IHC 1702 – EXAMINATION COUCH

- CRCA rectangular tube frame with CRCA sheet cabinets
- Retractable step stool
- Fixed upholstered top
- Three cabinets and three drawers.
- Finish Pre-Treated and Epoxy Powder Coated
- B.P. apparatus tray

DIMENSIONS

Length : 1830 mm
Width : 610 mm
Height : 810 mm



IHC 1703 - EXAMINATION COUCH

- CRCA rectangular tube frame with CRCA sheet cabinets
- One drawer and one cabinet.
- Finish Pre-Treated and Epoxy Powder Coated
- Overall dimension : 1840mm L X 570mm W X 860mm H

DIMENSIONS

Length : 1830 mm
Width : 610 mm
Height : 810 mm



IHC 1711 - EXAMINATION TABLE

- CRCA tubular construction.
- Backrest adjustable on ratchet.
- Finish pre-treated and epoxy powder coated.

DIMENSIONS

Length : 1830 mm
Width : 525 mm
Height : 810 mm



IHC 1720 - NORMAL EXAMINATION TABLE

- CRCA tubular construction.
- Finish pre-treated and epoxy powder coated.

DIMENSIONS

Length : 1830 mm
Width : 525 mm
Height : 810 mm



CORPORATE OFFICE

Indtech Healthcare Pvt. Ltd, No. 6/858-M,
2nd Floor, Suite No.893, Valamkottil Towers
Judgemukku, Kakkanad, Kochi - 682021
Ph: +91 7034 111 182



CALICUT OFFICE

Indtech Healthcare Pvt. Ltd, Office 1513,
5th Floor, HiLITE Business Park, Calicut- 673014
Ph: +91 7560 822 822



FACTORY

Plot No: B5, Industrial Development Plot,
Paloor Kotta, Palachode (P.O), Angadippuram (Via),
Malappuram-679 338, Kerala, Ph: +91 9495 351458



ISO CERTIFIED 13485:2016 & 9001:2015



Customer Care

+91 7594 822 822

www.indtechhc.com

info@indtechhc.com