



Patient Care Furniture



ISO CERTIFIED 13485:2016 & 9001:2015



Indtech expertise in manufacturing innovative, comfortable and easy-to-use healthcare furnitures with utmost safety for patients in mind. Working for frontline healthcare workers including doctors and nurses and patients who seeks care and comfort in their difficult times, Indtech never compromise quality.

Our products ranges from ICU advanced care and ICU support to Dialysis, Examination support, Ward care, General support, Transfer support and Mother care. The right talent in production applied upfront allows us to meet almost all of our client's requirements in time.

QUALITY POLICY

The quality policy of IndTech Health Care Pvt. Ltd., is to provide products that meet or exceed customer expectations, comply with all applicable regulations, continually improve in everything we do, and maintain the effectiveness of the quality management system.

WHAT MAKES US SPECIAL



CUSTOMISED PRODUCTS

We provide our customers with the products that meet their specific demands and requirements.



ISO 13485:2016 CERTIFIED

The products are accredited to ISO 13485:2016 & 9001: 2015 and CE certification for European Economic Area



INTERNATIONAL STANDARDS

Every time we make sure that the products we manufacture meets the international export standards.



PRODUCT MAINTENANCE

We don't quit once the product is handed over to our clients. We provide maintenance services for the products on site.



ON TIME DELIVERY

We make sure that the products we manufacture get delivered to our clients on time.



STATE-OF-THE-ART PRODUCTION SYSTEM

We use the latest technologies to manufacture the products and ensure the utmost quality.

In the face of laborious working environment in hospitals, Indtech offers comfortable and easy to handle furnitures with utmost safety for patients. With varieties of different designs and grades from deluxe to normal, we meet requirements for each and every hospitals. Check out our products to know more.





LEARN MORE ABOUT INDTECH PATIENT CARE FURNITURES.

ICU Advanced Care
Dialysis
Ward Care
Transfer Support
Mother Care

ICU Support
Examination Support
General Support



IHC 1505 - MAYO'S TROLLEY(SS)

- Height adjustable by knob system
- S.S Tubular frame
- Four 100 mm Dia Castors with High grade synthetic body



IHC1506 - INSTRUMENT TROLLEY(SS)

DIMENSIONS (A)

Length : 690 mm

Width : 460 mm

Height : 810 mm

DIMENSIONS (B)

Length : 760 mm

Width : 510 mm

Height : 810 mm



IHC1507 DRESSING TROLLEY(SS)

DIMENSIONS

Length : 760 mm

Width : 510 mm

Height : 810 mm



IHC1841 - LINEN TROLLEY WITH SHELVES(SS)



IHC1840 - LINEN TROLLEY(SS)



IHC 1807 - FOOD TABLE (ADJUSTABLE BY GAS SPRING)

- Rectangular frame work Mounted on 50mm Dia casters
- ABS Top adjustable by pneumatic gas spring
- Finish Pre-Treated and Epoxy Powder Coated



IHC 1808 - FOOD TABLE (ADJUSTABLE BY GEAR)

- Rectangular frame work Mounted on 50mm Dia casters
- ABS Top adjustable by Gear mechanism
- Finish Pre-Treated and Epoxy Powder Coated



IHC 1809 - FOOD TABLE (ADJUSTABLE BY KNOB)

- Rectangular frame work Mounted on 50mm Dia casters
- Wooden Top adjustable by knob mechanism
- Finish Pre-Treated and Epoxy Powder Coated



IHC1851 - BED SIDE SCREEN



IHC 1801 - BED SIDE LOCKER (DELUXE)



IHC 1802 - BED SIDE LOCKER



IHC 1825 - REVOLVING STOOL



IHC 1825C - REVOLVING STOOL



IHC1835 - SALINE STAND



IHC1836 - HEAVY BASE SALINE STAND



IHC 1829 - FOOT STEP MS



CORPORATE OFFICE

Indtech Healthcare Pvt. Ltd, No. 6/858-M,
2nd Floor, Suite No.893, Valamkottil Towers
Judgemukku, Kakkanad, Kochi - 682021
Ph: +91 7034 111 182



CALICUT OFFICE

Indtech Healthcare Pvt. Ltd, Office 1513,
5th Floor, HiLITE Business Park, Calicut- 673014
Ph: +91 7560 822 822



FACTORY

Plot No: B5, Industrial Development Plot,
Paloor Kotta, Palachode (P.O), Angadippuram (Via),
Malappuram-679 338, Kerala, Ph: +91 9495 351458



ISO CERTIFIED 13485:2016 & 9001:2015



Customer Care

+91 7594 822 822

www.indtechhc.com

info@indtechhc.com