



PRODUCT CATALOGUE

PROTECT YOUR PATIENTS,
PRESERVE YOUR INSTRUMENTS.



Hospital Autoclave

- Sterilization at standard 121°C & 134°C cycles
- PLC based automation, with 6 pre programmed cycles and additional user defined cycles
- Robust SS construction for door, chamber, jacket and steam generator
- On-line printer for continuous monitoring by operator
- Automatic sliding doors sealed with the help of a silicone door gasket
- The gasket is pressurized by compressed-air and retracted by special vacuum mechanism
- State of the art SCADA (Supervisory Control & Data Acquisition) software allows remote operation & monitoring of the machine
- Reports generated by the software in
 - 1) Detailed Data Format
 - 2) Single page summary
 - 3) Graph {Trend} can be printed Digital Data Back-up facility
- Audio Alarm & Visual Display of Error messages on Video Screen and on the computer Screen too
- Supervisory access codes and Password protection for prevention of unauthorized access

