

PRODUCT CATALOGUE

**PROTECT YOUR PATIENTS,
PRESERVE YOUR INSTRUMENTS.**



Tabletop Autoclave {Class B}

- Chamber and door fabricated out of SS 316L grade steel
- Power supply is 220 V Ac/50Hz
- Digital data back-up facility
- Control and Printing: Fully Automatic PLC control system provides on-line documentation through in-built printer

Pre-Vacuum phase efficiently removes air, resulting in the full penetration of steam into wrapped, porous and hollow loads. The post -vacuum, drying phase, ensures complete drying of porous loads and hollow instruments. This guarantees that even the most difficult loads will be sterilized completely

Safety:

- A cycle will not start unless the door is properly locked
- Door will not open after completion of the cycle until atmospheric pressure and safe temperature are achieved inside the chamber
- Automated double locking safety device prevents door from opening while chamber is pressurized.

